



Equipment Decontamination Form

Equipment Make/Model Information & UTEP Tag #: _____

This signed form must be affixed to all UTEP laboratory equipment before any maintenance, repairs, or re-location can occur. This includes laboratory equipment that needs is to be picked-up by surplus, or returned to the manufacturer.

Check completed tasks below (if applicable); Sign & date:

- ☐ Remove and clear the interior and surface of any bottles, material, liquid or equipment.
- ☐ If potentially biohazardous or biohazardous materials were used in the equipment: Decontaminate with Cavicide or freshly prepared 10% Bleach Solution followed by 70% alcohol wipe.
- ☐ Clean equipment interior and exterior with a damp cloth and mild detergent.
- ☐ Ensure equipment is free of any radioactive sources/contamination; verify decontamination by use of a certified Geiger Counter.
- ☐ Remove any UTEP logo stickers or "Hazard" stickers not initially placed by manufacturer.
- ☐ Other: _____

Signature: _____

Date: _____

Print Name: _____

Rm. # & Lab Phone: _____

Upon inspection, this equipment is deemed as being decontaminated by:

EH&S: _____

Date: _____

Print Name: _____

Rm. # & Office Phone: _____